

Fixation of Parasteatoda embryos for conventional whole-mount in situ hybridization (WISH)

by Hiroki Oda and Yasuko Akiyama-Oda (JT Biohistory Research Hall) Updated April 5, 2019

<https://www.brh2.jp>

- 1) Put eggs in a dry glass-made dish.
- 2) Add 1-2ml of 100% bleach and swirl eggs using a Pasteur pipet. Be careful not disrupting eggs. During this step, which we usually take for about 3 min, the surface of each egg should become transparent.
- 3) Add 1-2 ml of water and swirl eggs to let them sink down.
- 4) Remove the solution and wash the eggs with water several times.
- 5) Remove the water from the dish.
- 6) Add a part of the fixative and heptane, which should be prepared in a glass vial in advance. Eggs must be settled at the interface of the two phases.
- 7) Transfer the eggs to the glass vial containing 3ml of fixative and 3ml of heptane.
- 8) Shake the vial at 150 rpm at RT for 2 h or more.
- 9) Optional: continue to shake it at 4°C overnight (or for several hours).
- 10) Remove the fixative thoroughly.
- 11) Add 3ml of PEM and transfer the eggs to a glass-made dish.
- 12) Remove the PEM as much as possible and add PEM+0.1%tween 20 (PEMt).
- 13) Dissect out embryos using forceps.
- 14) Transfer the embryos to a 1.5ml Eppen tube.
- 15) Incubate the embryos in a graded series of MetOH/PEMt solutions (25% , 50%, 75%, and 95%) and finally in 100% MetOH.
- 16) Change the solution to 100% EtOH.
- 17) Store the embryos at -20°C until use.

Alternative method:

- 13) Incubate the eggs in a graded series of MetOH/PEM solutions (25%, 50%, 75%, and 95%) and finally in 100% MetOH for 10-30 min each.
- 14) Dissect embryos in MetOH using forceps. [The eggs can be stored at -20°C before dissection.]
- 15) Transfer the embryos to a 1.5ml Eppen tube.
- 16) Change the solution to 100% EtOH.
- 17) Store the embryos at -20°C until use.

2x PEM stock:		Fixative:		
200mM PIPES		2x PEM	1.5mL	
2mM EDTA		Water	1.05mL	
2mM MgSO ₄		37% Formaldehyde	0.45mL	
pH 6.9		Total	3.0mL	